**Artisan™  
Jenner Wright Giemsa Stain Kit****Code / Réf. / Code-Nr. AR308****ENGLISH****Intended use**

For In Vitro Diagnostic Use.

Jenner Wright Giemsa Stain Kit is intended for laboratory use to differentiate, by light microscopy, hematopoietic cells in bone marrow tissue samples.<sup>1-2</sup> Routinely processed samples (paraffin-embedded) **cut at 2-3 microns** may be used. The preferred fixative is neutral buffered formalin. The clinical interpretation of any staining, or its absence, should be complemented by morphological studies and proper controls, and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist. All reagents in the Jenner Wright Giemsa Stain Kit have been optimally prepared for use on Dako's Artisan Link Staining System and require no mixing or diluting.

**NOTE:** The use of harsh acid for decalcifying and/or excessive time in acid may cause inadequate staining results in bone marrow biopsy samples.

**Summary and explanation**

The Jenner Wright Giemsa Stain Kit is typically used to identify hematopoietic cells (i.e., eosinophils, polymorph nuclear leukocytes, etc.) in bone marrow tissue.

**Principle of procedure**

Using light microscopy, hematopoietic cells are visible in Hematoxylin and Eosin preparations, but appear homogenous with other cells. Greater differentiation is achieved with the application of this special stain.

The Jenner Wright Giemsa Stain Kit is optimized for use on the Artisan Link Staining System with or without the Artisan Clearing Solution. Wash steps follow all of the staining steps. Following staining, the air dried slides are removed from the instrument, cleared with xylene or xylene substitute and mounted with appropriate mounting media. Results are interpreted using a standard light microscope.

**Control Tissue:** Bone Marrow/Spleen

**Reagents provided**

The Jenner Wright Giemsa Stain Kit is composed of five ready-to-use dispenser packs, containing:

**Methanol****METHANOL**

Methanol

**Jenner Solution****JENNER SOLUTION**

Jenner Dye in methanol

**Wright Giemsa Solution****WRIGHT GIEMSA  
SOLUTION**

Wright/Giemsa in methanol

**Buffer Solution pH 6.75****BUFFER SOLUTION  
pH 6.75**

Phosphates in deionized water

**Acetic Acid 1%****ACETIC ACID 1%**

Acetic Acid in deionized water

Dako provides these reagents at optimal formulation for use according to the instructions provided for staining on prepared tissue sections. Any deviation from the recommended test procedures may invalidate declared expected results.

### Materials required, but not supplied

1. Microscope slides
2. Artisan Link Staining System
3. Wash Solution (code AR102)
4. Artisan Clearing Solution (code AR309) Optional

### Precautions

1. For professional users.
2. Minimize microbial contamination of reagents or increase in nonspecific staining may occur.
3. Reagent packs may contain residual solution after completing the prescribed number of tests.
4. Incubation times or temperatures other than those specified may give erroneous results. The user must validate any changes made to the factory-released procedure.
5. Store reagent packs in an upright position when not in use on the instrument.
6. Keep reagent packs out of direct sunlight and away from heat-generating sources.
7. As a general rule, persons under 18 years of age are not allowed to work with this product. Users must be carefully instructed in the proper work procedures, the dangerous properties of the product and the necessary safety instructions. Please refer to Safety Data Sheet (SDS) for additional information.
8. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
9. Unused solution should be disposed of according to local, State and Federal regulations.
10. Safety Data Sheet available for professional users on request.
11. For the vials shipped in hazmat bags, please check to see if there is any spilled liquid or broken reagent packs. If either of these conditions is present, do not open the hazmat bag. Return reagent pack to the kit box and report damage to customer service for a replacement kit.

**Danger**

#### **Methanol:** 60-100% Methanol

- H225 Highly flammable liquid and vapor.  
H301 + H311 + Toxic if swallowed, in contact with skin or if inhaled.  
H331  
H319 Causes serious eye irritation.  
H315 Causes skin irritation.  
H360 May damage fertility or the unborn child.  
H370 Causes damage to organs. (central nervous system (CNS))  
H335 May cause respiratory irritation.  
P201 Obtain special instructions before use.  
P202 Do not handle until all safety precautions have been read and understood.  
P281 Use personal protective equipment as required.  
P280 Wear protective gloves. Wear eye or face protection. Wear protective clothing.  
P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.  
P241 Use explosion-proof electrical, ventilating, lighting and all material-handling equipment.  
P242 Use only non-sparking tools.  
P243 Take precautionary measures against static discharge.  
P233 Keep container tightly closed.  
P271 Use only outdoors or in a well-ventilated area.  
P260 Do not breathe vapor.  
P270 Do not eat, drink or smoke when using this product.  
P264 Wash hands thoroughly after handling.  
P307 + P311 IF exposed: Call a POISON CENTER or physician.  
P304 + P340 + IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Call a POISON CENTER or physician.  
P311  
P301 + P310 + IF SWALLOWED: Immediately call a POISON CENTER or physician. Rinse mouth.  
P330  
P303 + P361 + IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.  
P353  
P302 + P352 + IF ON SKIN: Wash with plenty of soap and water. Call a POISON CENTER or physician if you feel unwell. Take off  
P312 + P362-2 contaminated clothing.  
P332 + P313 If skin irritation occurs: Get medical attention.  
P305 + P351 + IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.  
P338 Continue rinsing.  
P337 + P313 If eye irritation persists: Get medical attention.  
P405 Store locked up.  
P403 Store in a well-ventilated place.  
P235 Keep cool.  
P501 Dispose of contents and container in accordance with all local, regional, national and international regulations.

**Danger**

#### **Jenner Solution:** 60-100% Methanol

- H225 Highly flammable liquid and vapor.

H301 + H311 + H331 Toxic if swallowed, in contact with skin or if inhaled.  
H319 Causes serious eye irritation.  
H315 Causes skin irritation.  
H360 May damage fertility or the unborn child.  
H370 Causes damage to organs. (central nervous system (CNS))  
H335 May cause respiratory irritation.  
P201 Obtain special instructions before use.  
P202 Do not handle until all safety precautions have been read and understood.  
P281 Use personal protective equipment as required.  
P280 Wear protective gloves. Wear eye or face protection. Wear protective clothing.  
P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.  
P241 Use explosion-proof electrical, ventilating, lighting and all material-handling equipment.  
P242 Use only non-sparking tools.  
P243 Take precautionary measures against static discharge.  
P233 Keep container tightly closed.  
P271 Use only outdoors or in a well-ventilated area.  
P260 Do not breathe vapor.  
P270 Do not eat, drink or smoke when using this product.  
P264 Wash hands thoroughly after handling.  
P307 + P311 IF exposed: Call a POISON CENTER or physician.  
P304 + P340 + P311 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Call a POISON CENTER or physician.  
P301 + P310 + P330 IF SWALLOWED: Immediately call a POISON CENTER or physician. Rinse mouth.  
P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.  
P302 + P352 + P312 + P362-2 IF ON SKIN: Wash with plenty of soap and water. Call a POISON CENTER or physician if you feel unwell. Take off contaminated clothing.  
P332 + P313 If skin irritation occurs: Get medical attention.  
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P337 + P313 If eye irritation persists: Get medical attention.  
P405 Store locked up.  
P403 Store in a well-ventilated place.  
P235 Keep cool.  
P501 Dispose of contents and container in accordance with all local, regional, national and international regulations.

**Danger**

**Wright Giemsa Solution:** 60-100% Methanol, 0.1-1% Giemsa's stain, 0.1-1% Stains, biological, Wright's

H225 Highly flammable liquid and vapor.  
H301 + H311 + H331 Toxic if swallowed, in contact with skin or if inhaled.  
H319 Causes serious eye irritation.  
H315 Causes skin irritation.  
H360 May damage fertility or the unborn child.  
H370 Causes damage to organs. (central nervous system (CNS))  
H335 May cause respiratory irritation.  
P201 Obtain special instructions before use.  
P202 Do not handle until all safety precautions have been read and understood.  
P281 Use personal protective equipment as required.  
P280 Wear protective gloves. Wear eye or face protection. Wear protective clothing.  
P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.  
P241 Use explosion-proof electrical, ventilating, lighting and all material-handling equipment.  
P242 Use only non-sparking tools.  
P243 Take precautionary measures against static discharge.  
P233 Keep container tightly closed.  
P271 Use only outdoors or in a well-ventilated area.  
P260 Do not breathe vapor.  
P270 Do not eat, drink or smoke when using this product.  
P264 Wash hands thoroughly after handling.  
P307 + P311 IF exposed: Call a POISON CENTER or physician.  
P304 + P340 + P311 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Call a POISON CENTER or physician.  
P301 + P310 + P330 IF SWALLOWED: Immediately call a POISON CENTER or physician. Rinse mouth.  
P303 + P361 + P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water or shower.  
P302 + P352 + P312 + P362-2 IF ON SKIN: Wash with plenty of soap and water. Call a POISON CENTER or physician if you feel unwell. Take off contaminated clothing.  
P332 + P313 If skin irritation occurs: Get medical attention.  
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do.

P338	Continue rinsing.
P337 + P313	If eye irritation persists: Get medical attention.
P405	Store locked up.
P403	Store in a well-ventilated place.
P235	Keep cool.
P501	Dispose of contents and container in accordance with all local, regional, national and international regulations.

### Storage

Store at 2-8 °C. Do not use after expiration date. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, appropriate tissue must be run as controls for such verification. If unexpected staining is observed which cannot be explained by variations in laboratory procedures or improper storage, contact Dako Technical Support.

### Instrument procedure / Staining procedure

#### Preliminary Comments

Please refer to the *Artisan Link Basic User Guide* for detailed instructions.

#### Priming Reagent Packs

Before beginning staining, reagent packs should be primed and checked for leakage. Check the prime of each reagent pack when a new kit is used and before the first staining run of the day.

1. Over a paper towel, if present, remove the red or black shipping cap from the dispense assembly of the reagent pack by gently pulling down on the cap.
2. To prime a reagent pack when it is first received, hold the reagent pack with the dispenser assembly pointed downward and lightly tap on the side of the reagent pack 2–3 times to release any air bubbles that may be trapped in the dispenser assembly.
3. Fully compress and release the plunger 3–4 times in rapid succession, taking care to absorb any expelled fluid into a paper towel, absorbent gauze, or similar material, until a steady stream is obtained.
4. Avoid dispensing excess amounts of reagent during priming. Excessive priming will result in reagent loss, which may cause the reagent pack to run out of fluid before the stated number of tests is performed.
5. Before each use, reagent packs may be primed 1–2 times, checking for a steady stream.
6. Holding the reagent pack in one hand, slide the mounting guide on the back of the reagent pack into the slots on the reagent carousel until the pack is completely seated.

#### Procedure

Refer to the *Protocol* section in the Dako Link Software for the recommended detail procedure.

After staining is complete, slides will be dry. Dip slides in 95% Alcohol, 100% Alcohol, xylene or xylene substitute and permanently mount.

### Staining Interpretation

Nuclei:.....**Blue**  
 Eosinophils:.....**Bright Pink**  
 Leukocytes:.....**Shades of pink, gray, or blue depending on cell type and development**

### Troubleshooting

Refer to the Maintenance and Troubleshooting section of your *Artisan Link Basic User Guide* or contact our Technical Support.

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Dako North America, Inc.
Dako Denmark A/S